

Cimbria HEID GmbH

Heid Werkstrasse 4

2000 Stockerau - Austria

Telefon: (43) 2266-699-0

Telefon: (43) 2266-699101

Telefax: (43) 2266-65590

e-mail : heid@cimbria.at or wbo@cimbria.at

Internet: www.cimbria.com or www.cimbriakaack.com



Laboratory Equipment

Sample huller type Limprimita with husk separator



ELECTRIC SAMPLE HULLER WITH HUSK SEPARATOR

In special execution with reduction gear motor of 0,25 kW, combined with hull separator (winnow) with 30 watts single phase motor, huller and husk separator fitted on a common steel frame, the complete set readily assembled.

DECORTIQUEUR ELECTRIQUE D' ECHANTILLONS AVEC SEPARATEUR DE COSSES

En exécution spéciale avec moteur réducteur de 0,25 kW, combiné avec séparateur de cosSES (vanneur) avec moteur de 30 watts monophasé, décortiqueur et séparateur de cosSES assemblé sur un cadre commun en acier, l'unité complète prête assemblée.